

# HeartCheck™ PALM

Handheld ECG Devices  
Powered by CardioComm Solutions, Inc.

## The HeartCheck™ Palm Cardiac Event Monitor


Capture Lead I or II heart rhythm recordings.

Connect to a smart device with the GEMS™ Mobile App via Bluetooth to import ECG recordings.

Transmit ECGs to a physician, clinic or the SMART monitoring ECG reading service or generate your own ECG record in PDF format to show to your physician.

Available OTC the HeartCheck™ Palm allows a user to record a medical grade 30 second single lead ECG. Coupled with the GEMS™ Mobile App a user can import their recording, enter specific symptom, activity and diary entry information related to the recording and send for a review by a physician, clinic or the SMART monitoring ECG reading service.

Capable of fitting into multiple clinic and telemedicine workflows, the HeartCheck™ Palm with the GEMS™ Mobile App is an ideal tool for remote patient monitoring in today's ever expanding telemedicine market.



For more information please contact  
CardioComm Solutions, Inc. at

[sales@cardiocommsolutions.com](mailto:sales@cardiocommsolutions.com) or +1-416-977-9425 ex. 1



## Key Features

FDA Cleared, CE mark

Bluetooth Protocols:  
V4.0, Classic and Low Energy

Battery:  
Rechargeable lithium battery

Lead: I and II

Calibration voltage: 1mV ±5%

Standard sensitivity: 10 mm/mV ±5%

Scanning speed: 25 mm/s ±5%

Sampling rate: 250 hz

HR measurement range:  
30bpm ~ 300bpm,  
error ±1 bpm or 1%

Degree of waterproof: IP22

Display: 1.77" color TFT-LCD

Size: 100 mm (L) \* 45 mm (w) \* 15  
mm (H)

Weight: 60g



**CardioComm  
Solutions Inc.**

259 Yorkland Road, Suite 200  
North York, ON Canada